

## BT-100 Electric Breast Pump

# HI-bebe plus

**Anytime Anywhere quietly for mother and baby**



- **Bright and clear 3" LCD display**
- **Memory customized expression sequence**
- **Vacuum control in 10 levels**
- **Auto switch off after 30 min**
- **Sufficient vacuum for dual expression**
- **Backflow prevention**
- **Long time use rechargeable battery (150min)**
- **Convenient and graceful design**
- **3 operating mode (express, massage, memory)**
- **Operating time display**
- **Silent**

- **Intended use and indications of breast pumps**

Breast pumps are intended to be used by lactating women for expressing breast milk in home environments. Expressing milk by breast pump is indicated to alleviate the symptoms of engorgement. In case of symptoms of mastitis, breast pumps may support the healing process by removing breastmilk from the affected breast.

- **BT-100 can be used as following cases**

- ✓ Separated from baby because of work or travel
- ✓ To relieve sore and cracked nipples
- ✓ To bring out flat or retracted nipples
- ✓ When baby cannot suck mother's nipple because of severe nipple retraction
- ✓ To store the breastmilk because the milk expression amount is irregular

- **Strength**



## Specifications : BT-100 Breast Pump

Functional Characteristics			
Type		Pumping Pressure	
Dual Electric Breast Pump		Vacuum	-20 ~ -270mmHg
Display		Function	
Type	3" Mono LCD	Mode	Expression 10levels Massage 10levels
Pressure Level	Numeric & 10levels	Auto Shut Off	30min
Operating Time	Numeric Display(min:sec)	Memory	1 ~ 60levels
Battery Level	3levels	Back Flow Prevention	
Power			
Adaptor	Input : AC 100 ~ 240V (50/60Hz) Output : DC 12V/1.5A	Rechargeable Battery	7.4V Li-ion Polymer 900mA Operating Time : 150min Charging Time : 90min
Consumption	13.4W		
Standard Configuration			
Operation manual	1ea	Funnel & Bottle Set	1ea (Funnel : 24mm)
Power Adaptor	1ea		
Options			
Carrying Bag		Funnel & Bottle Set (for Dual)	
Funnel	21mm / 27mm	Cigar Cable	
USB Charging Cable	USB to Adaptor		
Warranty			
Main Unit	1year		
Physical Characteristics			
Dimension		Weight	
Main Unit	95(W) x 49(D) x 132(H)mm	Main Unit	260g
Packing	233(W) x 156(D) x 181(H)mm	Packing	0.8Kg
Carton Box (10ea)	805(W) x 255(D) x 385(H)mm	Carton Box (10ea)	10Kg
Environmental Conditions			
Operating Temperature	10 ~ 40°C (50 ~ 104°F)	Storage Temperature	-20 ~ 60°C (-4 ~ 140°F)
Operating Humidity	5 ~ 85% non-condensing	Storage Humidity	0 ~ 95% non-condensing
Certificates			
Certificates	KFDA, CE	Waterproof	IP21
Standard	IEC60601-1, IEC60601-1-2		

## Components

